


Maintenance

The ProTouch for Mighty Mouse protector is designed to capture finger oils, food particles and dirt that can accumulate over time. To remove built-up grime, we recommend that you remove, and then thoroughly wash the protector in warm water infused with mild dishwashing detergent. Regular washing may also help prevent the spread of germs. Dry thoroughly before reapplying the ProTouch to your Mighty Mouse.

▲ **Dry thoroughly before reapplying to the Mighty Mouse.**



Send us your feedback

We'd love to hear from you. Customer comments and suggestions are important elements in the development of future iSkin products.

Send your feedback to: feedback@iSkin.com

For FAQs regarding this product, visit us online at iSkin.com



iSkin and ProTouch are trademarks of iSkin Inc., registered in Canada, the U.S. and other countries. MICROBAN is a registered trademark of Microban Products Company. Apple and the Apple logo are trademarks of Apple Inc. Other company and product names may be trademarks of their respective owners.

1 Year Limited Warranty

iSkin Inc. warrants to the original purchaser that the product will be free from defects in materials and/or workmanship for the period of 1 Year. The warranty period commences on the date of purchase from the authorized reseller/retailer. Your original purchase invoice, showing the date of purchase of the product from the reseller/retailer is your proof of the date of purchase. This proof of purchase is mandatory for any warranty issues. iSkin Inc. at its option will repair or replace, at no charge, any part that iSkin Inc. determines to be defective.

This warranty shall not apply if the product has been damaged by abuse, misuse, negligence, accident, modification, tampering or by any other causes unrelated to defective materials and/or workmanship including acts of God. iSkin's warranty only covers product(s) made by iSkin and does not cover third party product(s) sold by other companies. iSkin acknowledges that with normal use of products made by iSkin, the purchased product(s) will provide satisfactory performance as specified by product performance literature.

Repair or replacement as provided under this warranty is iSkin's exclusive right. ANY APPLICABLE IMPLIED WARRANTIES, INCLUDING WARRANTIES OR MERCHANTABILITY AND FITNESS FOR A PARTICULAR PURPOSE ARE HEREBY LIMITED TO 365 DAYS FROM THE DATE OF PURCHASE AND ARE SUBJECT TO CONDITIONS SET FORTH IN THIS LIMITED WARRANTY. IN NO EVENT SHALL iSkin Inc. BE RESPONSIBLE FOR INCIDENTAL DAMAGES RESULTING FROM THE BREACH OF ANY EXPRESS OR IMPLIED WARRANTIES. This warranty gives you specific legal rights, and you may also have other legal rights, which vary from state to state, province to province, valid in the United States and Canada. Any shipping or other cost such as customs incurred costs, etc. from outside North America is the customers' responsibility.

Exclusions and Limitations

This Limited Warranty applies only to hardware product manufactured by iSkin Inc. that can be identified by the "iSkin" trademark, trade name, or logo affixed to it. This Limited Warranty does not apply to any non-iSkin products. This will include third party products packaged or sold with iSkin, any third party manufacturer may supply their own product warranty.

All inquiries with regards to the product warranty may be directed to: www.iSkin.com/warranty

iSkin

ProTouch

User's Guide




A The ProTouch Advantage

Made from premium silicone and featuring built-in antimicrobial properties to inhibit the growth of odor and stain causing bacteria on the product, the iSkin ProTouch for Mighty Mouse creates a barrier to keep dust, dirt and food from reaching your mouse. It also helps prevent surface wear, which can extend the usage life of your mouse.


iSkin ProTouch's premium silicone is long-lasting and can be easily washed and re-used as necessary. Regular washing helps maintain a clean surface, and can also help promote a healthier environment.

1 Prepare your Apple Mighty Mouse




Before using the iSkin ProTouch Mighty Mouse protector, unplug and wipe your mouse thoroughly with a damp, lint-free cloth to remove any surface oil and dirt that may have built up. Doing so will also ensure a better fit for the protector.

2 Overlaying the ProTouch protector




1 Place the ProTouch mouse protector over the Mighty Mouse. Start at the top end, ensuring that the mouse cable is slotted into the groove on the ProTouch.



2 Continue covering the body of the Mighty Mouse, ensuring that it is properly fitted, and that the scroll ball is aligned with the opening on the protector.

Overlaying the ProTouch protector contd...



3 Lift the outside rim of the ProTouch over the edges of the mouse. Align it with the contours of the mouse until it is properly fitted. Ensure that the scroll ball is fully exposed.

? About Microban® antimicrobial protection

Research studies have shown that the average keyboard can contain substantial amounts of bacteria. The iSkin ProTouch keyboard protector features a premium silicone material with a built-in antimicrobial product protection to inhibit the growth of odor and stain causing bacteria on the product - making it the first ever keyboard protector to offer this level of protection.



Notice to user:

iSkin Inc. shall not be responsible for any damage, including consequential damage caused by defective goods, personal injuries, property damage or consequential economic loss. iSkin does not guarantee suitability of its product for your application. Furthermore, iSkin Inc., shall not be responsible for any damage, loss of data, service, or otherwise caused to the computer or Mighty Mouse when used with this product.